



ATMOS^{Finland}**group**

CSSD Solutions

Full system to label, scan & track surgical equipment



ATMOS^{Finland}group

Atmos Group is a Finnish technology company manufacturing laser-marked labels for surgical instruments and accessories. We also develop custom-made instrument tracking systems for CSSDs to track items using our labels.

Brochure Content:

Pages 3 - 5: AtmosDot Labels

Pages 6 - 8: AtmosTrack Software

Page 9: Contact information





AtmosDot

Find and sort instruments faster with **color coded** labels.

Versatile Compatibility

Suitable for most materials.

Effortless Readability

Faster identifying and sorting.

Long-Lasting Durability

Last over 200 cycles.

Removable Design

Effortless removal when necessary.

Customizable Color Options

Choose your colors.

Easy Application

Apply and go!



General

AtmosDot labels are manufactured to order. The most common and recognized sizes are round, with diameters: 4.5 mm, 5.0mm and 6.5 mm.

Rectangular labels are also available, for long and narrow instruments.

Upon request, 2D matrix code (ECC200) can be laser marked on each label.
All matrix codes meet GS1 standards.

Each color has clearly contrasting markings for superior machine vision and human eye readability.

AtmosDot labels can be integrated in to any tracking system in use.





Customization

We understand that every central sterile services department (CSSD) has its own set of requirements and preferences when it comes to instrument tracking. That's why we offer fully customizable tracking labels, allowing you to tailor the shape, size, and markings to suit your specific needs.

The chart shows the PMS color of the labels and the outcome of the laser marked text / symbols / code.

186	ABC	286	ABC	724	ABC	black	ABC
347	ABC	PMS VIOLET	ABC	709	ABC	cool gray 8	ABC
124	ABC	543	ABC	528	ABC	1205	ABC
138	ABC	570	ABC	492	ABC	397	ABC
299	ABC	602	ABC	light gray	ABC	204	ABC



AtmosTrack

Validate wash loads, secure packaging & ensure sterilization.

Individual item tracking

Unique datamatrix codes for items.

Easy database creation

User friendly & fast item creation.

Easy integration

Cloud based SAAS solution.

Intuitive user interface


Just click & scan.

Cycle monitoring

Washing and autoclave validation & reports.

Enchanted communication

Web interface for OR.




General

While assembling sets in CSSD, labeled instruments are scanned into the AtmosTrack system, providing real-time location and OR data. This immediate data capture allows staff to instantly verify the presence of all necessary components within a set, minimizing errors and delays.

Furthermore, the system simultaneously updates the status of each instrument, preventing potential shortages and optimizing the flow of sterile supplies to the operating rooms.

The software aids correct container packing and helps locate replacement instruments. Built-in images show assembled trays and guide multi-piece instrument assembly.





Always have the right instruments, at the right place, at the right time

Key benefits

- Increased operational efficiency and reduced redundancy
- Improved compliance with health and safety regulations
- Data-driven insights for continuous process improvement
- Reduction of overload
- Maintenance and repair scheduling
- Streamlined workflow
- Enhanced communication between operating rooms and sterile processing



Contact information



Chief Executive Officer | Markku Korhonen
mk@atmosgroup.fi

Over 20 years of experience in the development and consulting of instrument sterilization solutions.



Chief Operating Officer | Elton Korhonen
elton@atmosgroup.fi

Served as CEO and project coordinator for multiple companies.

Office

Atmos Group Finland LTD
Sullituntie 5, 14700 Hauho, Finland
+358 40 681 7575
mk@atmosgroup.fi

Reg.nbr: FI34396007

